


ISO CERTIFIED • LOW PROFILE FLOOR SCALE SUMMIT 3000 SERIES

SUMMIT THE ECONOMICAL CHOICE OFFERING EXCEPTIONAL PERFORMANCE.



The **DIGI SUMMIT 3000 SERIES** provides exceptional performance in light industrial and legal-for-trade applications. This rugged, heavy duty industrial design provides outstanding durability in medical applications leaving the competition far behind. The deck will not bend or bow eliminating stress failures and increasing accuracy.

While value priced, the Summit 3000 also employs a rugged design to deliver efficient installation and operations. The convenient top access junction box is housed in a durable ABS for access ease.

The scale has fully adjustable, rigid polyurethane rubber feet to accommodate weighing on uneven surfaces.

STANDARD FEATURES

- Four NTEP - certified alloy steel shear beam load cells
- Top access junction box
- Rugged, heavy duty diamond safety treadplate steel deck
- Four adjustable rigid polyurethane rubber feet 3" (76.2 mm) diameter
- Single trim summing board housed in ABS enclosure prevents failure due to moisture and intrusions
- One Threaded eyebolt hole located in the center of the weldment for easy lifting
- 20 ft. of EL 147 load cell cable
- Conduit for additional load cell cable protection

INDICATORS - IQ355

- Flush or wall mount displays
- 8" light Emitting Diode (LED) display
- Stainless steel washdown enclosure
- 5 button operation.
- Digital paper print out (Patients weight and date)




OPTIONS/ACCESSORIES

Pit Frame - for in-floor scales



Access Ramps - for easy access for above ground applications



SUMMIT 3000 SERIES

Product #	Size	Capacity
Sum3x3	3x3	1000lb (1000 kg)
Sum4x4	4x4	1000lb (1000 kg)
Sum5x5	5x5	1000lb (5000kg)

Total Quality Medical, Inc.

Your Choice for Quality Dialysis Equipment & Supplies