



CODE 45010189

Giselle® Maintenance Cartridges-kit WW

Contains 4 clean cartridges package.

Package size: 100x100x100mm (WxDxH)

Cartridge Weight: 45gx4

Cartridge size: 40mmxØ75mm

High density polyethylene cartridge, sealed with protective film and child safety cap.

Ingredient: each cartridge contains 25g of citric acid monohydrate 100% for performing cleaning of Giselle®.

Cartridge content is identified by RFID label.



For more information please contact us at:

Total Quality Medical, Inc.

128 York Street, Kennebunk, ME 04043 • tqmedical@aol.com

Tel 207-985-5590 • 800-300-8776 • www.TQMedical.com