

TQM

Total Quality Medical



CODE 45000990

Premium cart

Close lockable cart, specifically design to house and fix Giselle® on the top and store capacity for tanks and other accessories such as salt cartridges, spray bottles, printer in the lower part.

The cart is equipped with shelves, key lockable doors, 4 non-marking grey wheels and 2 independent breaks.

Thanks to the cart Giselle® and its accessories can be easily moved without disassembly and reassembly.

*Materials: main structure - anodized aluminium; wheels - rubber caster
External dimensions: 570x585x920mm (WxDxH)
Wheel diameter: 100mm
Load-bearing capacity per shelf: 20kg*



For more information please contact us at:

Total Quality Medical, Inc.

128 York Street, Kennebunk, ME 04043 • tqmedical@aol.com

Tel 207-985-5590 • 800-300-8776 • www.TQMedical.com