

# TQM

Total Quality Medical



CODE 45000982

## Soleva<sup>®</sup> Spray Bottles Kit\_USA

Contains 2 blue Soleva<sup>®</sup> spray bottles + 2 green Soleva<sup>®</sup> spray bottles.

Blue Soleva<sup>®</sup> bottle:

- Designed for 0.1% free available chlorine concentration
- Size: 80x96x290mm (WxDxH)
- Weight (empty): 120g
- Blue connector and nebulizer for producing and applying Soleva<sup>®</sup> solution
- Blue label includes the solution's instruction for use

Green Soleva<sup>®</sup> bottle:

- Designed for 0.6% free available chlorine concentration
- Size: 64x95x300mm (WxDxH)
- Weight (empty): 120g
- Green connector and nebulizer for producing and applying Soleva<sup>®</sup> solution
- Green label includes the solution's instruction for use

Bottles are made of high-density polyethylene. Available space for application of manual or printed label with solution's production date and time.



For more information please contact us at:

**Total Quality Medical, Inc.**

128 York Street, Kennebunk, ME 04043 • [tqmedical@aol.com](mailto:tqmedical@aol.com)

Tel 207-985-5590 • 800-300-8776 • [www.TQMedical.com](http://www.TQMedical.com)