



CODE 45000989
Standard cart

Light open dismantable cart, specifically design to house and fix Giselle® on the top shelf and storage the two tap and waste water tanks on the lower shelf.

Thanks to the cart the equipment can be easily moved without disassembly and reassembly.

The cart is equipped with 2 shelves, 4 non-marking grey wheels and 2 independent breaks.

Materials: main structure - anodized aluminium; wheels - rubber caster
External dimensions: 540x540x1085mm (WxDxH)
Wheel diameter: 75mm
Load-bearing capacity per shelf: 20kg