

TQE Plus Endotoxin Filters

US FDA 510k Cleared

TQM
Total Quality Medical, Inc.

TQE Plus

Cartridge filter for Dialysis
Central Water System



Assists in preparing water to meet the
ANSI/AAMI/ISO 11663:2015 standard.

The TQE Plus Endotoxin Filters are designed for use in Dialysis Central Water Systems. The 5nm pore size offers the tightest control of endotoxins, bacteria, and viruses available today.

Product Features

- Available in 10", 20" and 30" length configurations
- Retains endotoxins and other contaminants greater than 5nm
- Expected life up to 6 months
- High flow rate capacity

Specifications


Specifications	TQE Plus
Compatible Disinfection	Chemical
Max Differential Pressure	30 psi (2 bar)
Filter Membrane ¹	Medisulfone®
Membrane Material	Polysulfone
Pore Size	5 nm
Bacteria Retention	> 10 ¹¹ (B. diminuta)
Virus Retention	> 10 ⁸ (PhiX-174)
Endotoxin Retention	> 10 ⁵ EU/ml
Sizes	10"/20"/30" x 2.5" D
Membrane surface area	2.6 m ² per 10" segment
Connection	222 O-Ring
Replacement ²	up to 6 months

1. Medisulfone® is a registered trademark of Medica S.p.A.

2. Up to 6 months with weekly chemical disinfection



10" TQE Plus - Clean Water Flow



Part Number	Description
TQE0.005-10P3E	TQE Plus Endotoxin 10" Filter Flat 222 Head
TQE0.005-20P3E	TQE Plus Endotoxin 20" Filter Flat 222 Head
TQE0.005-30P3E	TQE Plus Endotoxin 30" Filter Flat 222 Head

Distributed by:



Total Quality Medical, Inc.

128 York Street, Kennebunk, ME 04043
 Tel 207-985-5590 Toll Free 800-300-8776 Fax 207-985-5837
 sales@tqmedical.com www.TQMedical.com