



# MEDCOR™

"Medical Chair Specialist"

## Standard Features:

- 350 lb Patient Weight Capacity
- Two Folding Side Tables
- Four 5" swivel locking thermal plastic casters
- Trendelenburg position activated from either side of chair
- Built-in lumbar support
- Swing-away side arms for ease of cleaning & lateral transfer of patients
- All steel frame



## Optional Features:

- 450 / 550 lb. patient weight capacity
- Power Recline
- Heat and Massage
- Heat Only
- Headrest
- Headrest Cover
- Arm Pillows
- 3 Inch Leg Extensions

## Warranty

Chair frame	10 years
Reclining mechanism	2 years
Upholstery	2 years



### Model 180

### Model 350 R

#### SPECIFICATIONS

Seat Height	21"
Seat Depth	20"
Reclined Depth	76"
Seat Width	22"
Overall Width with side tables out	49.5"
Overall Width with side tables down	31.5"
Weight	160 lbs.
Recommended Weight Capacity	350 lbs.

### Model 180 XL

### Model 450 XL-R

#### SPECIFICATIONS

Seat Height	21"
Seat Depth	20"
Reclined Depth	76"
Seat Width	26"
Overall Width with side tables out	53.5"
Overall Width with side tables down	35.5"
Weight	180 lbs.
Recommended Weight Capacity	450 lbs.

For more information contact:  
Becky Adkins or Amanda Berrier

Phone: 888-579-1050

E-Fax: 1-877-938-7307

amanda.medcor@gmail.com

### Model 550

### Model 550 XXL-R

#### SPECIFICATIONS

Seat Height	21"
Seat Depth	20"
Reclined Depth	76"
Seat Width	30"
Overall Width with side tables out	57.5"
Overall Width with side tables down	39.5"
Weight	200 lbs.
Recommended Weight Capacity	550 lbs.



180 Lay Flat

Trendelenburg

